



INSTALLATION SHEET

Can-Dan Wall Mount Change Table
Model # 19-88-62
Model # 19-88-72

Installation Requirements:

The following information on the wall construction **MUST** be filled out and returned to our office..

Specification:

Electric control height adjustable from 22" to 37" (from floor)
Mattress size available - 62" x 24" and 72" x 24"
Mattress has 2 body safety belts
Three (3) wire baskets are located under the table top for storage.

Location for Table:

Length of open wall space for the Change Table _____

Electrical outlet (110VAC) available close by? Yes No

NOTE: Electrical outlet must be minimum of 18" off centerline of table location.

The table will extend 32" from the wall. Is there any interference with a door. Yes No

Type of Wall Construction:

1) Wood Studs:

Size of Studs: _____" x _____" Spacing of Studs _____"

Covering Board Type _____ Thickness _____"

2) Steel Studs _____" x _____" Spacing of Studs _____"

Covering Board Type _____ Thickness _____"

Important Note: In most cases the Wall Mounted Change Table can not be mounted to a standard wall stud construction. The wall **MUST** be reinforced to manufacturer specifications.

Do you need a copy of the Wall Preparation Instruction Sheet? Yes No

3) Concrete Blocks

Size of Blocks _____" x _____" Height of Block Wall _____"

Thickness of Block _____"

4) Other Wall Construction: Please specify:

Optional Floor Frame:

The option Floor Frame allows the Wall Mounted Change Table to sit directly on the floor.

There is no connection to the wall.

If there is any concern with the strength of the wall or re-construction or reinforcement of the wall, is not feasible, then the Floor Mount should be considered.

The Change Table will sit an extra 4" out from the wall.

Is the Floor Frame the best solution for you?

Yes No

What type of floor? _____ Concrete? Yes No

Company: _____ Date: _____

Contact Name: _____ Phone # _____