

CUSTOM WALL MOUNT ELECTRIC HEIGHT ADJUSTABLE CHANGE TABLE
Part Number 20-00-62 (standard)
Part number 20-00-62-CU (custom length available from 39" to 62")



Measurements: 62" x 24" Table Size (standard). Custom length available from 39" to 62".
Mattress: 2" thick.
Height Adjustment: 20" to 36" (without mattress.)
Minimum Length for space: 65" (standard). For custom lengths, add 3" additional clearance to table length.
Minimum Width Space Requirement: 33" (Folded 10").
Lifting Capacity: 250 lbs (evenly distributed).
Optional: 40" Side rail, 1-side only, outside table side, position offset on table.
Side Rail flips down for folding.
Weight of unit: 310 lbs.

Product Group: Change Tables

Description: The Wall Mount Changing Table (folding) is a lifting device and ETL Certified. Electrically powered hand control connected to a 24V actuator. The table is lowered for safe transfer for both user and caregiver. Individual can be secured using the 2 integrated safety belts. Table is raised to caregiver's ideal working height. Height adjustable with wired hand control. Optional Side Rail is released to swing away to facilitate transfer with a locking pin. Suitable for use in schools, clinics, private homes and other caregiving situations for children.

Materials: Powder-coated steel frame (table & side rail).
18 oz. vinyl mattress, mildew resistant, latex-free lead-free, phthalate free and fire retardant to NFPA 701 USA and CAN/ULC S109 standards.
The front cover is made of acrylic for easy cleaning.

Specification Sheet

Colour:	Side rail – Grey Frame & Tabletop – White Mattress – Grey
Adjustability:	Height adjustable. Optional fold down, swing away and removeable Side Rail.
Ergonomics:	High quality powered change tables. Provide a stable platform for changing individuals. Ability to be lowered to a suitable level for safe transfer of the user. Optional Side Rail offset on table to support and safely position user on table. Facilitates ease of access by caregiver.
Functionality:	Solid surface with a 2" pad. Flat surface facilitates assisted or independent transfers which are comfortable and safe. Hand control can be fastened to front cover. Removable mattress fastened with Velcro. Caregivers use the hand control to lower the change table to a suitable height to safety transfer individual to the padded tabletop.
Actuator:	LA31 Careline, (Linak)
Control Box:	CB9 Careline, (Linak)
Hand Control:	HB4 (Linak)
Motor:	Powered by 115VAC to 24 VDC power supply. Requires a GFI receptacle.
Cleaning:	Must be wiped down with a soft rag immediately after cleaning to prevent discolouration of the mattress. If not, the surface coating may deteriorate.
Warranty:	Can-Dan Rehatec Ltd. as manufacturer warrants from the DATE OF DELIVERY TO THE FIRST CONSUMER that each item described below, as NEW AND UNUSED, will be free from any defects in material and/or workmanship for a period of: <ul style="list-style-type: none"> • Four (4) years on change Table Frame (Manufactured by Can-Dan Rehatec Ltd.) • Two (2) year on Linak Actuator • One (1) year for accessories – accessories Include: <ul style="list-style-type: none"> · Hand Control · Tabletop · Mattress · Side rail · Safety Belts · Front Cover

Diagram with Dimensions (inches)

